



Your Partner for
**Reliable, Consistent,
Product Quality &
World Class Service**

DISPOSABLE MEDICAL DEVICE • PREVENTION
HEALTHCARE • SEXUAL WELLNESS • PLEASURE

Nulutex Sdn Bhd is a privately owned manufacturing company of latex condoms based in Kluang, Johor, Malaysia.

Established in 2004, Nulutex Sdn Bhd has produced and marketed quality latex condoms to over 40 countries around the world. Located close to the suppliers of our key raw material, i.e. natural rubber latex, we are able to offer our customers superior products at competitive price.

OUR MISSION

Nulutex Sdn. Bhd. is committed to manufacture and distribute excellent latex condoms and related personal care products that promote healthy lifestyle. It is also our aim to make latex condoms more popular and pleasurable.

OUR VISION

To be a global brand recognized for its quality disposable medical products, its contributions in creating a more caring society and its promotion of healthy lifestyle.

Managing Quality

In managing quality in our production, we emphasize and insist on using only the highest grade of raw materials. Our natural rubber latex*, which is the main raw material component in condom, is sourced from reputable suppliers. By using the specially formulated latex, we are able to produce remarkably thin condoms that are strong, reliable, and compliant with the required international safety and testing standards.

**Natural rubber latex complies with the requirements specified in the US Food and Drug Administration (FDA) Regulation CFR 21-884.5300*



OUR CORE VALUES



Reliable Quality



Competitive & Economic Cost



Customer Service



Product Safety



Flexibility & Customisation



Nulutex Sdn. Bhd.

Lot 1870, 4th Mile, Jalan Mengkibol, 86000 Kluang, Johor, Malaysia.

+607 777 1801

+607 777 1805

info@nulutex.com

www.nulutex.com



Latex Condoms

- See product range available below which can be packed in customised packaging and branding.



feel
feel+

House Brands

- Premium condoms for mass distribution complying to the WHO / UNFPA and ISO 4074 requirements.
- Properties:
Variants: Smooth / Studded / Ribbed
Widths / Sizes: 49mm / 53mm / 56mm
Other Option: Natural / Coloured / Flavoured



Nulatex Probe Covers

- Manufactured from highest quality natural rubber latex.
- Suitable for use on probes of Ultra Sound Imaging System. i.e. in trans-vaginal and trans-rectal imaging.
- Non-lubricated.
- Available non-sterile or sterile.
- Length = 200 mm
- Width = 53 mm
- Thickness = 0.066 to 0.075 mm



Intimate Lube

- Water based lubricant.
- Enhance pleasure with its silky feel on skin.
- Non-sticky and non-greasy.
- Safe to use with condoms.
- Customisable



Sachet 5g



Tube 50g / 100g



Bottle Pump 100g



Bottle Cap 100g

PRODUCTS	COLOUR	ABBREVIATION	WIDTH	LENGTH*	THICKNESS	SURFACE	DIAGRAM	
NORMAL RANGE	Smooth	○ Natural	NS	49 mm±2mm	180 ± 5 mm	66 µm - 75 µm	Regular Smooth	
		○ Natural	NS	53 mm+2mm	185 ± 5 mm			
		● Yellow	YS					
		● Violet	VS					
		● Red	RS					
		● Light Red	LRS					
		● Orange	OS					
		● Green	GS					
		● Pink	PS					
		● Blue	BLS					
		● Black	BS					
		● Brown	BRW					
	Extra Thick	○ Natural	NS			56 mm±2mm	195 ± 5 mm	90 µm - 110 µm
		● Light Red	LRS(ES)	53 mm+2mm	185 ± 5 mm			
		● Pink	PS(ES)					
Ultra Thin	○ Natural	NT	53 mm±2mm	185 ± 5 mm	56 µm - 65 µm	ULTRATHIN		
	● Light Red	LRT						
	● Pink	PT						
Probe Cover	○ Natural	NT	56 mm±2mm	195 ± 5 mm				
	○ Natural	NPC	53 mm+2mm	200 ± 5 mm	66 µm - 75 µm	Smooth		
TEXTURED & SHAPED RANGE	Studded	○ Natural	NST	49 mm±2mm	180 ± 5 mm	70 µm - 80 µm	Regular Textured	
		○ Natural	NST	53 mm±2mm	185 ± 5 mm			
		● Violet	VST					
		● Yellow	YST					
		● Pink	PST					
		○ Natural	NST					
	Ribbed	○ Natural	NR			53mm+ mm	185 ± 5 mm	70 µm - 80 µm
		● Pink	PR					
	3 in 1	○ Natural	N3C	53mm+2 mm	185 ± 5 mm	70 µm - 80 µm	Shaped Textured	
		● Pink	P3C					
		● Yellow	Y3C					
	Contour	○ Natural	NC	53mm+2mm	185 ± 5 mm	66 µm - 75 µm	Shaped Smooth	
		● Pink	PC					
	Flared	○ Natural	NF	53mm+2mm	195 ± 5 mm	66 µm - 75 µm	Flared	
	PREMIUM RANGE	Super Ultra Thin	○ Natural	SUT	53mm+2mm	185 ± 5 mm	45 µm - 55 µm	Extra Thin Smooth
Super Studded		○ Natural	NSST	53mm+2mm	185 ± 5 mm	90 µm - 110 µm	Extra Textured	