



**Nulatex Sdn Bhd** is a privately owned manufacturing company of latex condoms and related personal care products based in Kluang, Johor, Malaysia.

Established in 2004, Nulatex Sdn Bhd has produced and marketed quality latex condoms internationally. Located close to the suppliers of our key raw material, i.e. natural rubber latex, we are able to offer our customers products at competitive price.

### OUR MISSION

Nulatex is dedicated to manufacturing latex condoms and related personal care products that promote a healthy lifestyle. Our mission is to empower individuals with reliable, safe, and enjoyable options for sexual wellness, ensuring their well-being and satisfaction.

### OUR VISION

Our vision is to be the global leader in the manufacturing and distribution of disposable latex medical products, renowned for our commitment to innovation, quality, and customer satisfaction. We aspire to create a world where every individual has access to safe options for sexual wellness, fostering a culture of health and responsibility.

### Managing Quality

In managing quality in our production, we emphasize and insist on using only good raw materials. Our natural rubber latex\*, which is the main raw material component in condoms, is sourced from reputable suppliers. By using this good quality latex, we are able to produce remarkably thin condoms that are strong, reliable, and compliant with the required international safety and testing standards.

*\*Natural rubber latex complies with the requirements specified in the US Food and Drug Administration (FDA) Regulation CFR 21-884.5300*



### OUR CORE VALUES



Reliable Quality



Competitive & Economic Cost



Customer Service



Product Safety




Flexibility & Customisation



 Lot 1870, 4th Mile, Jalan Mengkibol, 86000 Kluang, Johor, Malaysia.

 +607 777 1801

 info@nulatex.com

 +607 777 1805

 [www.nulatex.com](http://www.nulatex.com)

 @feelbynulatex

 @feelbynulatex

 Nulatex Sdn Bhd



### OEM Latex Condoms

- See our available product range on the next page which can be packed in customised packaging and branding.

feel



### House Brands

- Latex condoms for mass distribution complying with the ISO 4074 requirements.
- Properties:  
 Variants: Smooth / Studded / Ribbed  
 Widths / Sizes: 49mm / 53mm / 56mm  
 Other Option: Natural / Coloured / Flavoured

feel

MDA Registration No.:  
GC93961343817



MDA Registration No.:  
GC38026230117



MDA Registration No.:  
GC7630620-44423



PRODUCTS	COLOUR	ABBREVIATION	FLAVOURS	CE APPROVED	FDA 510K	CFDA	WIDTH	LENGTH	THICKNESS	SURFACE	DIAGRAM
SMOOTH	○ NATURAL	NS	None	Yes	Yes	Yes	49 mm ± 2 mm	180 ± 5 mm	66 µm - 75 µm	Regular Smooth	
			Apple			Yes					
			Jasmin			Yes					
	● PINK	PS	None			Yes					
			Apple			Yes					
			Jasmin			Yes					
	○ NATURAL	NS	None	Yes	Yes	Yes	53 mm ± 2 mm	185 ± 5 mm			
			Apple			Yes					
			Jasmin			Yes					
			Vanilla	Yes	Yes						
	Strawberry	Yes									
	● YELLOW	YS	Banana	Yes	Yes						
	● VIOLET	VS	Grape	Yes							
	● RED	RS	Strawberry	Yes							
	● LIGHT RED	LRS		Yes							
	● ORANGE	OS	Orange	Yes	Yes						
● PINK	PS	None	Yes	Yes	Yes						
		Apple			Yes						
		Jasmin			Yes						
		Strawberry	Yes								
● BROWN	BRW	Chocolate	Yes								
○ NATURAL	NS	None	Yes			56 mm ± 2 mm	195 ± 5 mm				
			Yes			60 mm ± 2 mm	min 190 mm	50 µm - 80 µm			
EXTRA THICK	● PINK	NS (ES)	Yes			53 mm ± 2 mm	185 ± 5 mm	90 µm - 110 µm			
		PS (ES)	Strawberry	Yes							
ULTRA THIN	○ NATURAL	NT	None	Yes	Yes	Yes	53 mm ± 2 mm	185 ± 5 mm			
			Apple			Yes					
			Jasmin			Yes					
			Strawberry	Yes							
			Orange	Yes							
	Banana	Yes									
	● PINK	PT	None	Yes	Yes	Yes					
			Apple			Yes					
Jasmin					Yes						
Strawberry	Yes	Yes									
○ NATURAL	NT	None	Yes			56 mm ± 2 mm	195 ± 5 mm				
		Yes	Yes	Yes							
STUDDED	○ NATURAL	NST	Apple			Yes	53 mm ± 2 mm	185 ± 5 mm			
			Jasmin			Yes					
			None	Yes	Yes	Yes					
	● PINK	PST	Apple			Yes					
			Jasmin			Yes					
			Strawberry	Yes							
RIBBED	○ NATURAL	NR	None	Yes	Yes	Yes	53 mm ± 2 mm	185 ± 5 mm			
			Apple			Yes					
			Jasmin			Yes					
	● PINK	PR	None	Yes	Yes	Yes					
			Apple			Yes					
			Jasmin			Yes					
			Strawberry	Yes							